




# GOLDEN4 - diameters of lines - Edelrid - ( III )

MIDDLE OF THE CANNOPY



- Dyneema LIROS**  
 0.6 = LIROS DC 060  
 0.7 = LIROS DC 100  
 0.8 = LIROS DC 120
- COUSIN Vectran:**  
 0.6 = 12100
- COUSIN Dyneema:**  
 0.95 = COUSIN 0.95 mm - red  
 1.11 = COUSIN 1.1 mm - red
- Dyneema Edelrid:**  
 1.0 = 7850-080-018 (red)  
 1.0 = 7850-080-041 (yellow)  
 1.1 = 7850-100-041 (yellow)  
 1.3 = 7850-130-018 (red)  
 1.3 = 7850-130-041 (yellow)

- Aramid Edelrid:**  
 2.1 = 7343-280-018 (red)  
 2.1 = 7343-280-005 (bleu)  
 2.1 = 7343-280-041 (yellow)  
 1.8 = 7343-230-018 (red)  
 1.8 = 7343-230-005 (bleu)  
 1.8 = 7343-230-041 (yellow)  
 1.5 = 7343-190-018 (red)  
 1.5 = 7343-190-005 (bleu)  
 1.5 = 7343-190-041 (yellow)  
 1.3 = 7343-140-006 (green)

**BREAKS R**  
 Dyneema 2.0